

EXPERIENCE HIGH PERFORMANCE


PRECISION SPINE
SLIMPLICITY[®]HP
ANTERIOR CERVICAL PLATING SYSTEM

Dual thread cortical cancellous screws for
ROBUST SCREW PURCHASE

Individual level plate lordosis for
**OPTIMAL FIT & MINIMAL
PLATE BENDING**

Anatomically placed
bending zones to
EASE PLATE PLACEMENT

Single step blocker for
FAST & SIMPLE LOCKING



Fixed & variable screws to
MATCH ALL CERVICAL PROTOCOLS

Low Profile 2mm Plate to
**MINIMIZE TISSUE DISRUPTION &
POTENTIAL POST-OP DISCOMFORT**

Large Graft Windows for
**INCREASED INTRAOPERATIVE
VISUALIZATION
AND POSITIONING**

PRECISION SPINE[®]

Discover the Difference



SYSTEM TRAYS

P/N	Description
65-BK-0101	Plates & Screws + Standard Instruments
65-BK-0102	Plates & Screws + Large Handle Instruments

CERVICAL PLATES*

1 LEVEL PLATES

P/N	Description
65-CA-1010	1 Level - 10mm
65-CA-1012	1 Level - 12mm
65-CA-1014	1 Level - 14mm
65-CA-1016	1 Level - 16mm
65-CA-1018	1 Level - 18mm
65-CA-1020	1 Level - 20mm
65-CA-1022	1 Level - 22mm
65-CA-1024	1 Level - 24mm

2 LEVEL PLATES

P/N	Description
65-CA-2024	2 Level - 24mm
65-CA-2026	2 Level - 26mm
65-CA-2028	2 Level - 28mm
65-CA-2030	2 Level - 30mm
65-CA-2032	2 Level - 32mm
65-CA-2034	2 Level - 34mm
65-CA-2036	2 Level - 36mm
65-CA-2038	2 Level - 38mm
65-CA-2040	2 Level - 40mm
65-CA-2042	2 Level - 42mm
65-CA-2044	2 Level - 44mm

3 LEVEL PLATES

P/N	Description
65-CA-3043	3 Level - 43mm
65-CA-3046	3 Level - 46mm
65-CA-3049	3 Level - 49mm
65-CA-3052	3 Level - 52mm
65-CA-3055	3 Level - 55mm
65-CA-3058	3 Level - 58mm
65-CA-3061	3 Level - 61mm

4 LEVEL PLATES

P/N	Description
65-CA-4060	4 Level - 60mm
65-CA-4064	4 Level - 64mm
65-CA-4068	4 Level - 68mm
65-CA-4072	4 Level - 72mm
65-CA-4076	4 Level - 76mm
65-CA-4080	4 Level - 80mm

DISPOSABLES

FIXATION PINS

P/N	Description
65-FS-2010	Fixation Pin, Locking Mechanism
65-FS-2510	Fixation Pin, Screw Hole



PRECISION SPINE

SLIMPLICITY[®] HP

ANTERIOR CERVICAL PLATING SYSTEM

BONES SCREWS

Ø 4.0mm Self-Drilling FIXED SCREWS

P/N	Description
65-DF-4012	4.0mm x 12mm SD Fix Screw
65-DF-4014	4.0mm x 14mm SD Fix Screw
65-DF-4016	4.0mm x 16mm SD Fix Screw
65-DF-4018	4.0mm x 18mm SD Fix Screw

Ø 4.5mm Self-Tapping FIXED SCREWS






P/N	Description
65-TF-4512	4.5mm x 12mm ST Fix Screw
65-TF-4514	4.5mm x 14mm ST Fix Screw
65-TF-4516	4.5mm x 16mm ST Fix Screw
65-TF-4518	4.5mm x 18mm ST Fix Screw

Ø 4.0mm Self-Drilling VARIABLE SCREWS

P/N	Description
65-DV-4012	4.0mm x 12mm SD Var Screw
65-DV-4014	4.0mm x 14mm SD Var Screw
65-DV-4016	4.0mm x 16mm SD Var Screw
65-DV-4018	4.0mm x 18mm SD Var Screw

Ø 4.5mm Self-Tapping VARIABLE SCREWS

P/N	Description
65-TV-4512	4.5mm x 12mm ST Var Screw
65-TV-4514	4.5mm x 14mm ST Var Screw
65-TV-4516	4.5mm x 16mm ST Var Screw
65-TV-4518	4.5mm x 18mm ST Var Screw



*Plate Sizes measured center hole to center hole.
Add 8mm to plate length for end to end plate length

PRECISION SPINE[®]
www.precisionspineinc.com

