

EXPERIENCE OPTIMAL CONTROL

PRECISION SPINE
SURELOK™ MIS 3L
PERCUTANEOUS SCREW SYSTEM

Open tab design
ENHANCES VISUALIZATION AND EASES
ROD PLACEMENT

1.4mm low profile outer diameter
PROMOTES A SMALLER INCISION
& MINIMIZES MUSCLE DISRUPTION

Cannulated triple lead threads
INCREASES SPEED AND EFFICIENCY OF SCREW
DELIVERY

Proximal attachment feature
PROVIDES SECURE ATTACHMENT OF
STABILITY RINGS, SLEEVES AND
INTERNAL ROD REDUCTION
INSTRUMENTS

150mm integrated extension tab
ALLOWS 20MM OF CONTROLLED
ROD REDUCTION

30° of angulation in all planes
OFFERS INTRAOPERATIVE FLEXIBILITY

Self-starting aggressive screw tip
OPTIMIZES SCREW ENGAGEMENT




PRECISION SPINE®
Discover the Difference



SYSTEM TRAYS

P/N	Description
63-8600-CA-01	SureLOK MIS 3L Center Load Implants (6.5 & 7.5mm)
63-8700-CA-01	SureLOK MIS 3L Center Load Implants (5.5 & 8.5mm)
63-8300-CA-01	SureLOK MIS 3L Center Load Instruments Tray
63-8600-CA-02	SureLOK MIS 3L Side Load Implants (6.5 & 7.5mm)
63-8700-CA-02	SureLOK MIS 3L Side Load Implants (5.5 & 8.5mm)
63-8300-CA-02	SureLOK MIS 3L Side Load Instruments Tray



BONE SCREWS

Ø 5.5mm SCREWS

P/N	Description
63-CP-5535	MIS 3L Pedicle Screw - 5.5 x 35mm
63-CP-5540	MIS 3L Pedicle Screw - 5.5 x 40mm
63-CP-5545	MIS 3L Pedicle Screw - 5.5 x 45mm
63-CP-5550	MIS 3L Pedicle Screw - 5.5 x 50mm
63-CP-5555	MIS 3L Pedicle Screw - 5.5 x 55mm



Ø 6.5mm SCREWS

P/N	Description
63-CP-6535	MIS 3L Pedicle Screw - 6.5 x 35mm
63-CP-6540	MIS 3L Pedicle Screw - 6.5 x 40mm
63-CP-6545	MIS 3L Pedicle Screw - 6.5 x 45mm
63-CP-6550	MIS 3L Pedicle Screw - 6.5 x 50mm
63-CP-6555	MIS 3L Pedicle Screw - 6.5 x 55mm



Ø 7.5mm SCREWS

P/N	Description
63-CP-7535	MIS 3L Pedicle Screw - 7.5 x 35mm
63-CP-7540	MIS 3L Pedicle Screw - 7.5 x 40mm
63-CP-7545	MIS 3L Pedicle Screw - 7.5 x 45mm
63-CP-7550	MIS 3L Pedicle Screw - 7.5 x 50mm
63-CP-7555	MIS 3L Pedicle Screw - 7.5 x 55mm



Ø 8.5mm SCREWS

P/N	Description
63-CP-8535	MIS 3L Pedicle Screw - 8.5 x 35mm
63-CP-8540	MIS 3L Pedicle Screw - 8.5 x 40mm
63-CP-8545	MIS 3L Pedicle Screw - 8.5 x 45mm
63-CP-8550	MIS 3L Pedicle Screw - 8.5 x 50mm
63-CP-8555	MIS 3L Pedicle Screw - 8.5 x 55mm



RODS & LOCKING CAPS

LORDOTIC TI CENTER LOAD RODS

P/N	Description
63-LT-5035	Lordotic Ti CL Rod - 35mm
63-LT-5040	Lordotic Ti CL Rod - 40mm
63-LT-5045	Lordotic Ti CL Rod - 45mm
63-LT-5050	Lordotic Ti CL Rod - 50mm
63-LT-5055	Lordotic Ti CL Rod - 55mm
63-LT-5060	Lordotic Ti CL Rod - 60mm
63-LT-5065	Lordotic Ti CL Rod - 65mm
63-LT-5070	Lordotic Ti CL Rod - 70mm
63-LT-5075	Lordotic Ti CL Rod - 75mm
63-LT-5080	Lordotic Ti CL Rod - 80mm
63-LT-5090	Lordotic Ti CL Rod - 90mm
63-LT-5100	Lordotic Ti CL Rod - 100mm
63-LT-5110	Lordotic Ti CL Rod - 110mm



STRAIGHT TI CENTER LOAD RODS

P/N	Description
63-ST-5035	Straight Ti CL Rod - 35mm
63-ST-5040	Straight Ti CL Rod - 40mm
63-ST-5045	Straight Ti CL Rod - 45mm
63-ST-5050	Straight Ti CL Rod - 50mm
63-ST-5055	Straight Ti CL Rod - 55mm
63-ST-5060	Straight Ti CL Rod - 60mm
63-ST-5065	Straight Ti CL Rod - 65mm
63-ST-5070	Straight Ti CL Rod - 70mm
63-ST-5075	Straight Ti CL Rod - 75mm
63-ST-5080	Straight Ti CL Rod - 80mm
63-ST-5090	Straight Ti CL Rod - 90mm
63-ST-5100	Straight Ti CL Rod - 100mm
63-ST-5110	Straight Ti CL Rod - 110mm
63-ST-5120	Straight Ti CL Rod - 120mm
63-ST-5130	Straight Ti CL Rod - 130mm
63-ST-5140	Straight Ti CL Rod - 140mm
63-ST-5150	Straight Ti CL Rod - 150mm



NOTE

Side Load Rod part numbers:
Straight Ti Side Load Rods 48-ST-55XXX
Lordotic Ti Side Load Rods 48-CU-55XXX

LOCKING CAP

P/N	Description
SL1000	Cap Screw (Lock-Screw)



PRECISION SPINE®
www.precisionspineinc.com

